---

# PEEK

your space  
your air

---

## **HEPA-filtersysteem in BioSafetyLevel 3 en 4 omgevingen**

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**Peek b.v.  
Ruud Peek  
De Ring 13  
5261 LM VUGHT**

**16 april 2013**

**Referentienummer: RPeek/VOF1200012/V4**

Peek bv  
De Ring 13  
5261 LM Vught

T +31 (0)73 820 02 88  
F +31 (0)73 656 49 25  
E [info@peekbv.nl](mailto:info@peekbv.nl)

[www.peekbv.nl](http://www.peekbv.nl)  
KvK Tilburg 18030812  
ABN-AMRO 41.12.88.849

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<b>1. HEPA FILTERSYSTEEM IN BSL 3 EN 4 OMGEVINGEN.....</b>	<b>3</b>
<b>2. BSL 3 EN 4.....</b>	<b>3</b>
<b>3. UITGANGSDOCUMENT: NEN-EN 12 128:.....</b>	<b>7</b>
<b>4. NEN-EN 12 128: HFDST. 7.5/8.5: VENTILATION AND EXTRACTION .....</b>	<b>9</b>
<b>5. NEN-EN 12 128: HFDST. 7.7/8.7: SEALIBILITY .....</b>	<b>12</b>
<b>6. VEREISTE FILTERSYSTEEM CONFORM EN:12 128.....</b>	<b>14</b>
<b>7. UITBREIDING EISEN HEPA FILTERSYSTEMEN .....</b>	<b>15</b>
<b>8. TOTAALOVERZICHT HEPA FILTERSYSTEEM.....</b>	<b>20</b>

## 1. HEPA filtersysteem in BSL 3 en 4 omgevingen.

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- Biosafety level 1 (BSL 1)
- Biosafety level 2 (BSL 2)
- Biosafety level 3 (BSL 3)
- Biosafety level 4 (BSL 4)

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### Standards and Guidelines – EU

- EN 12 128 Laboratories for Research, Development and Analysis
- EN 1620 Biotechnology-Large Scale Process and Production
- EN 1822 HEPA and ULPA Filter

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### Standards and Guidelines – USA and Canada

- BMBL Biosafety in Microbiological and Biomedical Laboratories Centers of Disease Control and Prevention (CDC) & National Institutes of Health (NIH) 4th Edition April 1999
- Army Regulations and Pamphlet 385-69
- Canadian Laboratory Biosafety Guidelines Health Canada – 3rd Edition 2004

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**Table 1: Classification of Infectious Microorganisms by Risk Group**

<b>Risk Group Classification</b>	<b>NIH Guidelines for Research involving Recombinant DNA Molecules 2002<sup>2</sup></b>	<b>World Health Organization Laboratory Biosafety Manual 3<sup>rd</sup> Edition 2004<sup>1</sup></b>
Risk Group 1	Agents not associated with disease in healthy adult humans.	(No or low individual and community risk) A microorganism unlikely to cause human or animal disease.
Risk Group 2	Agents associated with human disease that is rarely serious and for which preventive or therapeutic interventions are often available.	(Moderate individual risk; low community risk) A pathogen that can cause human or animal disease but is unlikely to be a serious hazard to laboratory workers, the community, livestock or the environment. Laboratory exposures may cause serious infection, but effective treatment and preventive measures are available and the risk of spread of infection is limited.
Risk Group 3	Agents associated with serious or lethal human disease for which preventive or therapeutic interventions may be available (high individual risk but low community risk).	(High individual risk; low community risk) A pathogen that usually causes serious human or animal disease but does not ordinarily spread from one infected individual to another. Effective treatment and preventive measures are available.
Risk Group 4	Agents likely to cause serious or lethal human disease for which preventive or therapeutic interventions are not usually available (high individual risk and high community risk).	(High individual and community risk) A pathogen that usually causes serious human or animal disease and can be readily transmitted from one individual to another, directly or indirectly. Effective treatment and preventive measures are not usually available. <sup>3</sup>

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
**Table 2. Summary of Recommended Biosafety Levels for Infectious Agents**

BSL	Agents	Practices	Primary Barriers and Safety Equipment	Facilities (Secondary Barriers)
1	Not known to consistently cause diseases in healthy adults	Standard microbiological practices	<ul style="list-style-type: none"> <li>No primary barriers required.</li> <li>PPE: laboratory coats and gloves; eye, face protection, as needed</li> </ul>	Laboratory bench and sink required
2	<ul style="list-style-type: none"> <li>Agents associated with human disease</li> <li>Routes of transmission include percutaneous injury, ingestion, mucous membrane exposure</li> </ul>	BSL-1 practice plus: <ul style="list-style-type: none"> <li>Limited access</li> <li>Biohazard warning signs</li> <li>“Sharps” precautions</li> <li>Biosafety manual defining any needed waste decontamination or medical surveillance policies</li> </ul>	Primary barriers: <ul style="list-style-type: none"> <li>BSCs or other physical containment devices used for all manipulations of agents that cause splashes or aerosols of infectious materials</li> <li>PPE: Laboratory coats, gloves, face and eye protection, as needed</li> </ul>	BSL-1 plus: <ul style="list-style-type: none"> <li>Autoclave available</li> </ul>
3	Indigenous or exotic agents that may cause serious or potentially lethal disease through the inhalation route of exposure	BSL-2 practice plus: <ul style="list-style-type: none"> <li>Controlled access</li> <li>Decontamination of all waste</li> <li>Decontamination of laboratory clothing before laundering</li> </ul>	Primary barriers: <ul style="list-style-type: none"> <li>BSCs or other physical containment devices used for all open manipulations of agents</li> <li>PPE: Protective laboratory clothing, gloves, face, eye and respiratory protection, as needed</li> </ul>	BSL-2 plus: <ul style="list-style-type: none"> <li>Physical separation from access corridors</li> <li>Self-closing, double-door access</li> <li>Exhausted air not recirculated</li> <li>Negative airflow into laboratory</li> <li>Entry through airlock or anteroom</li> <li>Hand washing sink near laboratory exit</li> </ul>
4	<ul style="list-style-type: none"> <li>Dangerous/exotic agents which post high individual risk of aerosol-transmitted laboratory infections that are frequently fatal, for which there are no vaccines or treatments</li> <li>Agents with a close or identical antigenic relationship to an agent requiring BSL-4 until data are available to redesignate the level</li> <li>Related agents with unknown risk of transmission</li> </ul>	BSL-3 practices plus: <ul style="list-style-type: none"> <li>Clothing change before entering</li> <li>Shower on exit</li> <li>All material decontaminated on exit from facility</li> </ul>	Primary barriers: <ul style="list-style-type: none"> <li>All procedures conducted in Class III BSCs or Class I or II BSCs in combination with full-body, air-supplied, positive pressure suit</li> </ul>	BSL-3 plus: <ul style="list-style-type: none"> <li>Separate building or isolated zone</li> <li>Dedicated supply and exhaust, vacuum, and decontamination systems</li> <li>Other requirements outlined in the text</li> </ul>

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#### EN 12 128 Biotechnology – Laboratories for Research, Development and Analysis

Requirements	Physical contamination level			
	1	2	3	4
Containment level labelled	Yes <sup>2)</sup>	Yes	Yes	Yes
Hazard zones labelled with biohazard sign	-	Yes	Yes	Yes
Adequate space for each worker	Yes	Yes	Yes	Yes
Laboratory rooms separated by doors	No <sup>3)</sup>	Yes	Yes (lockable)	Yes (physically separated)
Entry to laboratory via an airlock	No	No	Optional	Yes
Provision of an observation window, or alternative	No	Yes	Yes	Yes
Surfaces impervious to water, easy to clean and resistant to cleaning agents etc.	Yes (bench)	Yes (bench)	Yes (bench, floor)	Yes (bench, floor, walls and ceiling)
Hand washing facilities	Yes	Yes	Yes	Yes
Provision of taps that can be operated without being touched by hand	Optional	Yes	Yes	Yes
Personnel shower	Optional	Optional	Optional	Yes (in airlock system)
Hand disinfection	Optional	Yes	Yes	Yes
Provision of storage for protective clothing within the laboratory or unit	Optional	Yes	Yes	Yes
<b>Ventilation</b> - Means of maintaining a negative pressure - Extract air to be filtered using HEPA filtration	No No	No No	Optional Yes	Yes Yes (via 2 HEPA filters)
<b>Input air to be filtered using HEPA filtration</b>	No	No	No	Yes
Alarm system fitted to detect unacceptable air pressure changes	No	No	Yes	Yes
Laboratory sealable to permit fumigation	No	Optional	Yes	Yes
Laboratory design to permit vector, e.g. insect and rodent, control	No	Optional	Yes	Yes
Laboratory to contain its own, dedicated, equipment	No	No	Yes	Yes
Provision of microbiological safety cabinet	Optional	Optional	Yes	Yes (class III)
Provision of emergency electrical supply	No	No	Optional	Yes
Means of outside communication, e.g. hands-free telephone or equivalent	No	No	Optional	Yes

<b>Provision of an autoclave</b> - accessible within laboratory suite - in laboratory, double-ended	No -	No -	Yes Optional	- Yes
<b>Waste disposal system</b> - documented - validated	No -	Yes -	- Yes	- Yes
<b>Effluent treatment</b>	No	Optional	Yes	Yes

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**BSL3:** All extract air shall be exhausted through a high efficiency particulate air (HEPA) filter. HEPA filters shall be sited so that they are accessible for testing and to allow their safe removal.

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**BSL4:** The exhaust air shall be filtered through two HEPA filters mounted in series before it leaves the laboratory unit.

*Note 2: The first HEPA filter should be sited within the laboratory*

*The input air shall be filtered through a single HEPA filter*

*HEPA filters shall be sited so that they are accessible for testing and to allow their safe removal.*

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(2) Filterpakking dichtheidstest

**Filter System SCF<sup>hightec</sup>**  
**Seal Test Device and Clamping Principle**

Together we can do it. **YIT**


**Filter System SCF<sup>hightec</sup>**  
**Seal Test Device and Clamping Principle**

1. Connection of leak test device (1) to fast acting coupling (2) at the door seat
2. Pumping air into the test groove by means of the hand pump (3) until the test pressure (2,000 Pa) is reached (readable at the pressure gauge (4))
3. Holding this pressure at a constant value by pressing the hand pump (3) smoothly
4. The readable volume (5) flow has to be lower than the admissible leakage rate

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**BSL 3 en 4:** *The laboratory and associated extract ventilation systems Shall be capable of being sealed for fumigation and shall allow effective removal or extraction of the fumigant.*

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
### (1) Gasdichte kleppen met meetril

		<p style="text-align: center;"><b>Gastight Dampers - Seal Test Device</b></p> <ol style="list-style-type: none"> <li>1. Connection of leak test device (1) to fast acting coupling (2) at the damper blade seat</li> <li>2. Pumping air into the test groove by means of the hand pump (3) until the test pressure (2,000 Pa) is reached (readable at the pressure gauge (4))</li> <li>3. Holding this pressure at a constant value by pressing the hand pump (3) smoothly</li> <li>4. The readable volume flow (5) has to be lower than the admissible leakage rate</li> </ol>  <p style="text-align: right; font-size: small;">Portable leak test device, type LTD</p> <p style="text-align: right; font-size: x-small;">Together we can do it. <b>YIT</b></p>
<p><b>Gasdichte klep met meetril</b></p>		

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



Auftragsnr. / manufacturer no.: G-8331			
<b>2. Tests Safe Change Filterhousing</b>			
<b>Type SCF TripleS 1x1 – H14</b>			
<b>2.1 leaktight testing filter housing</b>		<b>o.k.</b>	<b>failure</b>
2.1.1 HEPA- filter housing	see attachment 2.1.1	<input type="checkbox"/>	<input type="checkbox"/>
2.1.2 HEPA- filter housing, Triple S 1 <sup>st</sup> HEPA stage	see attachment 2.1.2	<input type="checkbox"/>	<input type="checkbox"/>
2.1.3 HEPA- filter housing, raw to clean gas side 1 <sup>st</sup> HEPA stage	see attachment 2.1.3	<input type="checkbox"/>	<input type="checkbox"/>
<b>2.2 seal testing filter seat and damper blade seat</b>		<b>o.k.</b>	<b>failure</b>
2.2.1 seal testing filter seat 1 <sup>st</sup> HEPA stage	see attachment 2.2.1	<input type="checkbox"/>	<input type="checkbox"/>
<b>2.3 clamping device HEPA stage</b>		<b>o.k.</b>	<b>failure</b>
check of clamping device by 3x release and clamp of filter element		<input type="checkbox"/>	<input type="checkbox"/>
<b>2.4 contamination free filter change</b>		<b>o.k.</b>	<b>failure</b>
function test of zero contamination filter exchange for HEPA-filter stage		<input type="checkbox"/>	<input type="checkbox"/>
<b>2.5 actuator of scanner system HEPA stage</b>		<b>o.k.</b>	<b>failure</b>
proper installation and wiring scanner is moving smoothly covering the whole filter area scanner actuator is stopped in correct position by limit switch „Top“ scanner actuator is stopped in correct position by limit switch „Bottom“		<input type="checkbox"/>	<input type="checkbox"/>
<b>2.6 completeness check</b>		<b>o.k.</b>	<b>failure</b>
check of completeness and marking according to contract specifications		<input type="checkbox"/>	<input type="checkbox"/>
<b>2.7 visual check</b>		<b>o.k.</b>	<b>failure</b>
visual control of welding and surfaces		<input type="checkbox"/>	<input type="checkbox"/>

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*Het copyright van dit document berust bij Technische Handelonderneming Peek B.V. Niets uit dit document mag worden veeveelvoudigd en/of openbaar gemaakt door middel van druk, fotokopie, microfilm, geluidsband, elektronisch of op welke andere wijze dan ook zonder voorafgaande schriftelijke toestemming van Technische Handelonderneming Peek B.V.*